

Product Datasheet

HIF1 alpha antibody (Biotin) (orb242410)

Description	HIF1 alpha antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Hypoxia-inducible factor 1-alpha protein (579-826AA)
Target	HIF1A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16665
Expiration Date	12 months from date of receipt.