

Product Datasheet

Hemoglobin antibody (HRP) (orb242420)

Description	Hemoglobin antibody (HRP)
Species/Host	Guinea pig
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Rabbit Hemoglobin (Native Protein)
Target	Hemoglobin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.