

Product Datasheet

amyloid A protein antibody (orb242421)

Description	amyloid A protein antibody
Species/Host	Rabbit
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Dog Serum amyloid A protein (19-129AA)
Target	SAA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19708
Expiration Date	12 months from date of receipt.