


Product Datasheet

Arabidopsis thaliana APX1 antibody (orb242429)

Description	Arabidopsis thaliana APX1 antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana L-ascorbate peroxidase 1, cytosolic protein (3-247AA)
Target	APX1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q05431
Expiration Date	12 months from date of receipt.



Western blot analysis of Arabidopsis tha...