

Product Datasheet

Baker's Chitin synthase 1 antibody (FITC) (orb242447)

Description	Baker's Chitin synthase 1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	FITC
Immunogen	Recombinant Baker's yeast Chitin synthase 1 protein (4-200AA)
Target	CHS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08004
Expiration Date	12 months from date of receipt.