

Product Datasheet

Canine distemper virus Fusion glycoprotein F0 antibody (FITC) (orb242463)

Description	Canine distemper virus Fusion glycoprotein F0 antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Canine distemper virus Fusion glycoprotein F0 protein (136-608AA)
Target	F
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P12569
Expiration Date	12 months from date of receipt.