


Product Datasheet

Ovalbumin antibody (Biotin) (orb242474)

Description	Ovalbumin antibody (Biotin)
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Chicken Ovalbumin (Native Protein)
Target	SERPINB14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01012
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): Ovalb...