

Product Datasheet

Albumin antibody (FITC) (orb242479)

Description	Albumin antibody (FITC)
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	FITC
Immunogen	Serum Albumin purified from Chicken serum (Native Protein)
Target	Wheat Gliadin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19121
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.