

Product Datasheet

Chlamydophila pneumoniae Major outer porin antibody (Biotin) (orb242486)

Description	Chlamydophila pneumoniae Major outer porin antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Chlamydophila pneumoniae Major outer membrane porin protein (24-389AA)
Target	ompA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27455
Expiration Date	12 months from date of receipt.