

Product Datasheet

Cyprinus carpio Serine protease inhibitor Kazal-type 2-like antibody (HRP) (orb242492)

Description	Cyprinus carpio Serine protease inhibitor Kazal-type 2-like antibody (HRP)
Species/Host	Rabbit
Reactivity	Fish
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Cyprinus carpio Serine protease inhibitor Kazal-type 2-like protein (18-77AA)
Target	Serine protease inhibitor Kazal-type 2-like protein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A0A059T4U4
Expiration Date	12 months from date of receipt.