

Product Datasheet

Epididymal secretory protein E1 antibody (FITC) (orb242499)

Description	Epididymal secretory protein E1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Canine
Conjugation	FITC
Immunogen	Recombinant Dog NPC intracellular cholesterol transporter 2 protein (22-149AA)
Target	NPC2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q28895
Expiration Date	12 months from date of receipt.