

Product Datasheet

k12 Pyruvate dehydrogenase antibody (HRP) (orb242508)

Description	k12 Pyruvate dehydrogenase antibody (HRP)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli Pyruvate dehydrogenase [ubiquinone] protein (1-572AA)
Target	poxB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07003
Expiration Date	12 months from date of receipt.