

Product Datasheet

Beta-lactamase TEM antibody (FITC) (orb242515)

Description	Beta-lactamase TEM antibody (FITC)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	FITC
Immunogen	Recombinant Escherichia coli Beta-lactamase TEM protein (24-286AA)
Target	bla
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62593
Expiration Date	12 months from date of receipt.