

Product Datasheet

NAD(P) transhydrogenase subunit beta antibody (HRP) (orb242532)

Description	NAD(P) transhydrogenase subunit beta antibody (HRP)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli NAD(P) transhydrogenase subunit beta protein (330-462AA)
Target	pntB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0AB67
Expiration Date	12 months from date of receipt.