

Product Datasheet

NAD(P)H-flavin reductase antibody (Biotin) (orb242533)

Description	NAD(P)H-flavin reductase antibody (Biotin)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli NAD(P)H-flavin reductase protein (2- 233AA)
Target	fre
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0AEN1
Expiration Date	12 months from date of receipt.