



Product Datasheet

GADD45B antibody (orb24255)

Description	GADD45B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Growth arrest and DNA damage-inducible protein GADD45 beta protein (1-160AA)
Target	GADD45B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O75293
Dilution Range	IHC-P:1:50-200, IF:1:50-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...



Immunofluorescence analysis of HepG2 cel...