


Product Datasheet

Hemoglobin antibody (FITC) (orb242555)

Description	Hemoglobin antibody (FITC)
Species/Host	Rabbit
Reactivity	Guinea pig
Conjugation	FITC
Immunogen	Guinea Pig Hemoglobin protein (Native Protein)
Target	Hemoglobin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): Guinea...