

Product Datasheet

Helicobacter pylori Orotidine 5'-phosphate decarboxylase antibody (Biotin) (orb242562)

Description	Helicobacter pylori Orotidine 5'-phosphate decarboxylase antibody (Biotin)
Species/Host	Rabbit
Reactivity	H. pylori
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Helicobacter pylori Orotidine 5'-phosphate decarboxylase protein (1-227AA)
Target	pyrF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P56155
Expiration Date	12 months from date of receipt.