



Product Datasheet

Alanine aminotransferase 1 antibody (orb242594)


Description	Alanine aminotransferase 1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Alanine aminotransferase 1 protein (2-496AA)
Target	GPT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24298
Dilution Range	WB:1:500-2000, IHC-P:1:20-500, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Rat heart tiss...



Western blot aanalysis of HepG2 whole ce...