

Product Datasheet

GCNT2 antibody (orb24263)

Description	GCNT2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human N-acetyllactosaminide beta-1,6-N-acetylglucosaminyl-transferase protein (201-402AA)
Target	GCNT2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NOV5
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.