

Product Datasheet

Breast cancer type 1 susceptibility protein antibody (orb242655)

Description	Breast cancer type 1 susceptibility protein antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Breast cancer type 1 susceptibility protein (24-255AA)
Target	BRCA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P38398
Expiration Date	12 months from date of receipt.