



Product Datasheet

BTB/POZ domain-containing protein 2 antibody (orb242659)


Description	BTB/POZ domain-containing protein 2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human BTB/POZ domain-containing protein 2 protein (321-525AA)
Target	BTBD2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BX70
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of 293T whole cel...



Immunofluorescent analysis of MCF-7 cell...