

Product Datasheet

CKI isoform epsilon antibody (Biotin) (orb242664)

Description	Rabbit polyclonal antibody against CKI isoform epsilon conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Casein kinase I isoform epsilon protein (1-416AA)
Target	CSNK1E
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49674
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.