



Product Datasheet

CKII subunit alpha' antibody (orb242667)


Description	CKII subunit alpha' antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Casein kinase II subunit alpha' protein (8-350AA)
Target	CSNK2A2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19784
Dilution Range	IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cel...



Immunofluorescence staining of PC-3 cell...



Immunohistochemistry of paraffin-embedde...