


Product Datasheet

Catechol O-methyltransferase antibody (orb242675)

Description	Catechol O-methyltransferase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Transmembrane O-methyltransferase protein (31-264AA)
Target	LRTOMT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WZ04
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of HepG2
cel...