

Product Datasheet

Collagen alpha-3(IV) chain antibody (Biotin) (orb242704)

Description	Collagen alpha-3(IV) chain antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Collagen alpha-3(IV) chain protein (1427-1668AA)
Target	COL4A3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01955
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.