

Product Datasheet

coronavirus OC43 Spike glycoprotein antibody (Biotin) (orb242712)

Description	coronavirus OC43 Spike glycoprotein antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human coronavirus OC43 Spike glycoprotein protein (15-344AA)
Target	S
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P36334
Expiration Date	12 months from date of receipt.