

Product Datasheet

DNA-directed RNA polymerase III subunit RPC10 antibody (orb242751)

Description	DNA-directed RNA polymerase III subunit RPC10 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human DNA-directed RNA polymerase III subunit RPC10 protein (1-108AA)
Target	POLR3K
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y2Y1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.