

Product Datasheet

DNA-directed RNA polymerase III subunit RPC3 antibody (HRP) (orb242758)

Description	DNA-directed RNA polymerase III subunit RPC3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human DNA-directed RNA polymerase III subunit RPC3 protein (1-534AA)
Target	POLR3C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BUI4
Expiration Date	12 months from date of receipt.