

Product Datasheet

Eukaryotic translation elongation factor 1 epsilon-1 antibody (Biotin) (orb242783)

Description	Eukaryotic translation elongation factor 1 epsilon-1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Eukaryotic translation elongation factor 1 epsilon-1 protein (2-174AA)
Target	EEF1E1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	043324
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.