



Product Datasheet

Eukaryotic translation initiation factor 3 subunit I antibody (orb242786)

Description	Eukaryotic translation initiation factor 3 subunit I antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Eukaryotic translation initiation factor 3 subunit I protein (1-323AA)
Target	EIF3I
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13347
Dilution Range	WB:1:2000-1:10000, IF:1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109(Lane 1), ...



Immunofluorescence analysis of HeLa cell...