

Product Datasheet

Eukaryotic translation initiation factor 5 antibody (HRP) (orb242793)

Description	Eukaryotic translation initiation factor 5 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Eukaryotic translation initiation factor 5 protein (1-431AA)
Target	EIF5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55010
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.