



Product Datasheet

Fibroblast growth factor 2 antibody (orb242806)


Description	Fibroblast growth factor 2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Fibroblast growth factor 2 protein (143-288AA)
Target	FGF2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09038
Dilution Range	WB:1:1000-1:5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse kidney & ...



Western blot analysis of HepG-2(Lane 1) ...



Immunohistochemical staining of Mouse lu...