



Product Datasheet

Fibronectin type 3 and ankyrin repeat domains protein 1 antibody (orb242815)

Description	Fibronectin type 3 and ankyrin repeat domains protein 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Fibronectin type 3 and ankyrin repeat domains protein 1 protein (1-345AA)
Target	FANK1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TC84
Dilution Range	IHC-P:1:20-1:200, IF:1:50-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human te...



Immunofluorescence analysis of MCF-7 cel...