

Product Datasheet

Gastric intrinsic factor antibody (Biotin) (orb242828)

Description	Rabbit polyclonal antibody against Gastric intrinsic factor conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Cobalamin binding intrinsic factor protein (19-417AA)
Target	CBLIF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27352
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.