

Product Datasheet

Guanylate cyclase activator 2B antibody (orb242854)

Description Guanylate cyclase activator 2B antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC

Immunogen Recombinant Human Guanylate cyclase activator 2B protein (37-112AA)

Target GUCA2B

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID [Q16661](#)

Dilution Range IHC-P:1:200-1:500, IF:1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.



Immunofluorescence analysis of HepG2 cells.



Immunohistochemical staining of human smooth muscle tissue.