



Product Datasheet

Haptoglobin antibody (orb242858)


Description	Haptoglobin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Haptoglobin protein (19-406AA)
Target	HP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00738
Dilution Range	WB:1:1000-1:5000, IHC-P:1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Mouse stomach ...



Immunohistochemistry of paraffin-embedde...