



Product Datasheet

Histone deacetylase 3 antibody (orb242882)


Description	Histone deacetylase 3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF, IHC, WB
Immunogen	Recombinant Human Histone deacetylase 3 protein (1-428AA)
Target	HDAC3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15379
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using Histone deacetylase 3 antibody.



Immunoprecipitation analysis of HeLa cells using Histone deacetylase 3 antibody.



Immunofluorescence analysis of HeLa cells using Histone deacetylase 3 antibody.