

Product Datasheet

HLA class II histocompatibility, DRB1-3 chain antibody (Biotin) (orb242898)

Description	HLA class II histocompatibility, DRB1-3 chain antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DRB1-3 chain protein (31-266AA)
Target	HLA-DRB1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01912
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.