

Product Datasheet

Inactive N-acetylated-alpha-linked acidic dipeptidase-like protein 2 antibody (Biotin) (orb242912)

Description	Inactive N-acetylated-alpha-linked acidic dipeptidase-like protein 2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Inactive N-acetylated-alpha-linked acidic dipeptidase-like protein 2 protein (143-795AA)
Target	NAALADL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q58DX5
Expiration Date	12 months from date of receipt.