

Product Datasheet

Interleukin-36 R antagonist antibody (HRP) (orb242938)

Description	Rabbit polyclonal antibody against Interleukin-36 R antagonist conjugated to HRP
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Interleukin-36 receptor antagonist protein (1-155AA)
Target	IL36RN
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UBH0
Expiration Date	12 months from date of receipt.