

Product Datasheet

Izumo sperm-egg fusion protein 4 antibody (Biotin) (orb242948)

Description	Izumo sperm-egg fusion protein 4 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Izumo sperm-egg fusion protein 4 protein (16-214AA)
Target	IZUMO4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q1ZYL8
Expiration Date	12 months from date of receipt.