

Product Datasheet

Keratin, type I cytoskeletal 14 antibody (FITC) (orb242953)

Description	Keratin, type I cytoskeletal 14 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Keratin, type I cytoskeletal 14 protein (1-472AA)
Target	KRT14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02533
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.