

Product Datasheet

immunoglobulin epsilon Fc R antibody (FITC) (orb242978)

Description	immunoglobulin epsilon Fc R antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Low affinity immunoglobulin epsilon Fc receptor protein (48-321AA)
Target	FCER2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06734
Expiration Date	12 months from date of receipt.