



Product Datasheet

LAG3 antibody (orb242984)


Description	LAG3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Lymphocyte activation gene 3 protein (29-450AA)
Target	LAG3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18627
Dilution Range	WB:1:500-5000, IHC-P:1:200-500, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...



Immunofluorescence analysis of HepG2 cel...



Immunohistochemical staining of human to...