

Product Datasheet

Membrane-associated tyrosine-and threonine-specific cdc2-inhibitory kinase antibody (orb243007)

Description	Membrane-associated tyrosine- and threonine-specific cdc2-inhibitory kinase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Membrane-associated tyrosine- and threonine-specific cdc2-inhibitory kinase protein (900-1101AA)
Target	PKMYT1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99640
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.