

Product Datasheet

N(4)-(beta-N-acetylglucosaminy)-L-asparaginase antibody (Biotin) (orb243041)

Description	N(4)-(beta-N-acetylglucosaminy)-L-asparaginase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human N(4)-(beta-N-acetylglucosaminy)-L-asparaginase protein (206-346AA)
Target	Genome polyprotein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20933
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.