

Product Datasheet

Nuclear nucleic acid-binding protein C1D antibody (FITC) (orb243065)

Description	Nuclear nucleic acid-binding protein C1D antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Nuclear nucleic acid-binding protein C1D protein (1-135AA)
Target	C1D
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13901
Expiration Date	12 months from date of receipt.