

Product Datasheet

Peptidyl-prolyl cis-trans isomerase FKBP3 antibody (FITC) (orb243088)

Description	Peptidyl-prolyl cis-trans isomerase FKBP3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP3 protein (2-224AA)
Target	FKBP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00688
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.