

Product Datasheet

Potassium voltage-gated channel subfamily A member 1 antibody (FITC) (orb243100)

Description	Potassium voltage-gated channel subfamily A member 1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily A member 1 protein (1-154AA)
Target	KCNA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q09470
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.